BPL_TEST2_Chemostat script with PyFMI

The key library PyFMI is installed.

After the installation a small application BPL_TEST2_Chemostat is loaded and run. You can continue with this example if you like.

```
!lsb_release -a # Actual VM Ubuntu version used by Google
```

```
No LSB modules are available.
Distributor ID: Ubuntu
Description: Ubuntu 22.04.4 LTS
Release: 22.04
Codename: jammy
```

%env PYTHONPATH=

```
→ env: PYTHONPATH=
```

!python --version

→ Python 3.11.11

```
!wget https://repo.anaconda.com/miniconda/Miniconda3-py311_24.11.1-0-Linux-x86_64.sh
!chmod +x Miniconda3-py311_24.11.1-0-Linux-x86_64.sh
!bash ./Miniconda3-py311_24.11.1-0-Linux-x86_64.sh -b -f -p /usr/local
import sys
sys.path.append('/usr/local/lib/python3.11/site-packages/')
```

!conda update -n base -c defaults conda --yes

Executing transaction: ...working... done

Platform: linux-64

installation finished.

Collecting package metadata (repodata.json): done

Solving environment: done

Package Plan

environment location: /usr/local

added / updated specs:
 - conda

The following packages will be downloaded:

package	build	
ca-certificates-2024.12.31 certifi-2025.1.31	 h06a4308_0 py311h06a4308_0	128 KB 163 KB
	Total:	291 KB

The following packages will be UPDATED:

Downloading and Extracting Packages:

```
BPL_TEST2_Chemostat_colab.ipynb - Colab
                                       | : 0% 0/1 [00:00<?, ?it/s]
| : 0% 0/1 [00:00<?, ?it/s]
| : 100% 1.0/1 [00:00<00:00, 10.06it/s]
    certifi-2025.1.31 | 163 KB
    ca-certificates-2024 | 128 KB
    certifi-2025.1.31 | 163 KB
    ca-certificates-2024 | 128 KB
                                         : 100% 1.0/1 [00:00<00:00, 8.30it/s]
    Preparing transaction: done
    Verifying transaction: done
    Executing transaction: done
!conda --version
!python --version
→ conda 24.11.1
     Python 3.11.11
!conda config --set channel_priority strict
    CustomValidationError: Parameter channel_priority = 'stric' declared in --set parameter is invalid.
     'stric' is not a valid ChannelPriority
    Valid choices for channel_priority: 'strict', 'flexible', 'disabled'
!conda \ install \ -c \ conda-forge \ pyfmi \ --yes \ \# \ Install \ the \ key \ package
₹
```

Preparing transaction: done Verifying transaction: done Executing transaction: done

Now specific installation run a simulation and notebook for that

Start with connecting to Github. Then upload the two files:

- FMU BPL_TEST2_Chemostat_linux_om_me.fmu
- Setup-file BPL_TEST2_Chemostat_explore.py

Filter out DepracationWarnings for 'np.float as alias' is needed — wish I could make filter more narrow import warnings warnings.filterwarnings("ignore")

%%bash

git clone https://github.com/janpeter19/BPL_TEST2_Chemostat

Cloning into 'BPL_TEST2_Chemostat'...

%cd BPL_TEST2_Chemostat

/content/BPL_TEST2_Chemostat

run -i BPL_TEST2_Chemostat_explore.py

→ Linux - run FMU pre-comiled OpenModelica

Model for bioreactor has been setup. Key commands:

- change of parameters and initial values - par()
- change initial values only - init()
- simulate and plot - simu() make a new plot - newplot()
- show()
- show plot from previous simulation
- disp() - display parameters and initial values from the last simulation - describe() - describe culture, broth, parameters, variables with values/units

Note that both disp() and describe() takes values from the last simulation


and the command process_diagram() brings up the main configuration

Brief information about a command by help(), eg help(simu) Key system information is listed with the command system_info()

```
%matplotlib inline
plt.rcParams['figure.figsize'] = [25/2.54, 20/2.54]
```

process_diagram()

No processDiagram.png file in the FMU, but try the file on disk.




describe('culture')


→ Simplified text book model – only substrate S and cell concentration X

```
newplot()
par(Y=0.50, qSmax=0.75, Ks=0.1)
                                              # Culture parameters
init(V_start=1.0, VX_start=1.0, VS_start=20) # Bioreactor startup
par(S_in=30, t0=0, F0=0, t1=10, F1=0.2)
                                              # Substrate feeding
simu(60)
```

Could not find cannot import name 'dopri5' from 'assimulo.lib' (/usr/local/lib/python3.11/site-packages/assimulo/lib/_i Could not find cannot import name 'rodas' from 'assimulo.lib' (/usr/local/lib/python3.11/site-packages/assimulo/lib/_in Could not find cannot import name 'odassl' from 'assimulo.lib' (/usr/local/lib/python3.11/site-packages/assimulo/lib/_i Could not find ODEPACK functions. Could not find RADAR5



```
# Check simu('cont')
newplot()
simu(20)
simu(40,'cont')
```



Start coding or generate with AI.