## → BPL\_TEST2\_Chemostat script with PyFMI

The key library PyFMI is installed.

After the installation a small application BPL\_TEST2\_Chemostat is loaded and run. You can continue with this example if you like.

```
!lsb_release -a # Actual VM Ubuntu version used by Google
    No LSB modules are available.
    Distributor ID: Ubuntu
    Description: Ubuntu 20.04.5 LTS
    Release:
                     20.04
    Codename:
                    focal
%env PYTHONPATH=
    env: PYTHONPATH=
!wget·https://repo.anaconda.com/miniconda/Miniconda3-py39_23.1.0-1-Linux-x86_64.sh
! \texttt{chmod} \cdot + \texttt{x} \cdot \texttt{Miniconda3-py39\_23.1.0-1-Linux-x86\_64.sh}
!bash../Miniconda3-py39_23.1.0-1-Linux-x86_64.sh.-b.-f.-p./usr/local.
import · sys ·
sys.path.append('_/usr/local/lib/python3.9/site-packages/')
    --2023-04-21 06:39:36-- https://repo.anaconda.com/miniconda/Miniconda3-py39_23.1.0-1-Linux-x86_64.sh
    Resolving repo.anaconda.com (repo.anaconda.com)... 104.16.131.3, 104.16.130.3, 2606:4700::6810:8203, ...
    Connecting to repo.anaconda.com (repo.anaconda.com) | 104.16.131.3 | :443... connected.
    HTTP request sent, awaiting response... 200 OK
    Length: 69888122 (67M) [application/x-sh]
    Saving to: 'Miniconda3-py39_23.1.0-1-Linux-x86_64.sh'
    Miniconda3-py39_23. 100%[===========] 66.65M
                                                               183MB/s
                                                                           in 0.4s
    2023-04-21 06:39:36 (183 MB/s) - 'Miniconda3-py39_23.1.0-1-Linux-x86_64.sh' saved [69888122/69888122]
    PREFIX=/usr/local
    Unpacking payload ...
    Installing base environment...
    Downloading and Extracting Packages
    Downloading and Extracting Packages
    Preparing transaction: done
    Executing transaction: done
    installation finished.
!conda update -n base -c defaults conda --yes
```

```
Preparing transaction: done Verifying transaction: done
Executing transaction: done \\
```

```
!conda --version
!python --version
```

conda 23.3.1 Python 3.9.16

!conda install -c conda-forge pyfmi --yes # Install the key package

```
Preparing transaction: done
Verifying transaction: done
Executing transaction: done
```

## Now specific installation run a simulation and notebook for that

Start with connecting to Github. Then upload the two files:

- FMU BPL\_TEST2\_Chemostat\_linux\_om\_me.fmu
- · Setup-file BPL\_TEST2\_Chemostat\_explore.py

```
# Filter out DepracationWarnings for 'np.float as alias' is needed - wish I could make filter more narrow
import warnings
warnings.filterwarnings("ignore")
%%hash
git clone https://github.com/janpeter19/BPL_TEST2_Chemostat
     Cloning into 'BPL_TEST2_Chemostat'...
%cd BPL_TEST2_Chemostat
     /content/BPL_TEST2_Chemostat
run -i BPL_TEST2_Chemostat_explore.py
    Linux - run FMU pre-comiled OpenModelica 1.21.0
     Model for bioreactor has been setup. Key commands:

    change of parameters and initial values
    change initial values only

     - par()
     - init()
     - simu() - simulate and plot

- newplot() - make a new plot
                   - show plot from previous simulation
      - show()
                    - display parameters and initial values from the last simulation
     - disp()
     - describe() - describe culture, broth, parameters, variables with values / units
    Note that both disp() and describe() takes values from the last simulation
     Brief information about a command by help(), eg help(simu)
    Key system information is listed with the command system_info()
%matplotlib inline
plt.rcParams['figure.figsize'] = [25/2.54, 20/2.54]
newplot()
par(Y=0.50, qSmax=0.75, Ks=0.1)
                                        # Culture parameters
init(V_0=1.0, VX_0=1.0, VS_0=20)
                                          # Bioreactor startup
par(S_in=30, t0=0, F0=0, t1=10, F1=0.2) # Substrate feeding
simu(60)
```

8

